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Reference

Tokyo Head Office

14F, JR Shinagawa East Bldg, 18-1, Konan 2-chome,  
Minato-ku, Tokyo 108-8230

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
※ Our membrane modules, the system with these modules installed and their related technical information may be subject to the regulation and approvals for specific applications by government.


For the exportation of these, you must have the necessary procedures like application of export license.

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Provide most suitable membrane · solution from laboratory to industrialization scale

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**食品分野** (Food field)

食品原料の精製・濃縮、飲料の清澄化 (Food ingredient refinery/concentration clarification of beverage)

詳細 ▶ p.12 (Details:p.12)

**浄水** (Purified water)

浄水処理 (UF中空糸) (Drinking water treatment (UF hollow fiber))

詳細 ▶ p.14 (Details:p.14)

**化学分野** (Field of chemistry)

染料・水溶性ポリマーの精製 (Refinery of dye and water soluble polymer)

詳細 ▶ p.13 (Details:p.13)

**海水** (Sea water)

海水無菌化 (Sea water sterilization)

詳細 ▶ p.14 (Details:p.14)

**医薬・バイオ** (Pharmaceutical/Bio)

酵素の菌体分離・濃縮、精製水製造 (Cell separation of enzyme concentration refinery water production)

詳細 ▶ p.13 (Details:p.13)

**排水** (Waste water treatment)

排水リサイクル、排水中の有価物回収 (Waste water recycle Recovery of valuable materials)

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MOLSEP® series of separation membrane products

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Polyethersulfone membrane of heat and chemical resisting and high separation performance.

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Polyacrylonitrile with less protein adsorption

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Long term water treatment reference

Our separation membrane "MOLSEP" series has begun with development of first domestic cellulose acetate made separation membrane for reverse osmosis in 1973.

As one of its application, our artificial dialysis water producing device "MOLSEP" which RO membrane applied has been popularized with spreading of dialysis treatment domestically.

And we have contributed greatly to reach the technique of membrane separation into the food industry filed through the juice concentration with no-heating and additive free process using our tubular type with cellulose acetate made RO membrane and success of commercial viability of fine sugar system taking advantage of tubular type membrane with heat-resisting type ultra filtration (Polyether sulfone made).

After that, we have succeed to develop the ultra filtration hollow fiber module with heat-resisting polyether sulfone membrane as well for refined water and ultra pure water production. Our ultra filtration hollow fiber membrane module with polyether sulfone made membrane developed by our original technique is highly evaluated at not only domestic but overseas due to its high rejection performance of endotoxin and its heat and chemical resisting character.

And also we had aimed to diversify of the hollow fiber membrane type and as a result they have been adopted to medical, pharmaceutical, and ultra pure water fields.

On the other hand, for the drinking water plant, we have succeed to develop the first domestic ultra filtration cellulose acetate made hollow fiber membrane and our large sized ultra filtration membrane module "FW50" was adopted at drinking water plant with membrane filtration method in 1995 as first domestic reference.

Our naturally derived ultra filtration cellulose acetate made hollow fiber membrane of hydrophilic character have been evaluated its high filtration performance and adopted widely.

Our separation membrane products "MOLSEP" series adopted widely, mainly domestic drinking water plant is widely used focuses of at the fields of food, medical, water treatment for long time in domestic and overseas.

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Cellulose acetate made hollow fiber membrane of high flux

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Abounding line ups in membrane materials, membrane module

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Superior heat and chemical resistance.

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Heat sterilization is possible due to its high heat resistance.

It is also excellent in chemical resistance, so it can be used in a wide variety of applications.

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Abundant in membrane molecular weight cut off

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
Abundant in membrane molecular weight cut off and hollow fiber diameter could be suitable with every applications.

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Superior lower protein adsorption

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
Hydrophilic polyacrylonitrile has less adsorption of protein character can maintain the stable filtration velocity.

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Sharp membrane separation performance

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Sharp molecular weight cut off curve

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\*1 It may be different depending on operating conditions. Especially, in case of constant high temperature use, please consult in advance.

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
Sustainable for long term of high flux operation.

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Cellulose acetate made hollow fiber membrane is unlikely to be fouling and can sustain the high flux for long term.

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
Behavior of flux by each membrane material

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By the routine chemical cleaning, membrane filtration performance will be recovered and realized the stable long term operation.

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Daily variation of trans membrane pressure difference of large membrane module

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Abundant experience for river and ground water treatment.

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Cellulose acetate made hollow fiber membrane module have been authorized as membrane module for drinking water production. With this module, you can sure about security/safe drinking water production. We have references more than 170 sites at domestic drinking water plants.

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スパイラル型 module  
Spiral type module

用途例詳細  
Application example Details

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MOLSEP® hollow fiber ultra filtration membrane module has high performance and prepared rich variations. Thanks to this, we can select the proper module size, membrane grade according to application and scale of the project and we can propose the most suitable modules based on water quality and your purpose.

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#### ■ Long-term stable operation possible

It is an internal pressure filtration structure. Cross flow filtration method enables long-term stable processing.

#### ■ High filtration performance

High pressure filtration performance can be maintained by back pressure washing periodically from the secondary side.

#### ■ Compactification of membrane treatment system.

Compact size because of hollow fiber membrane module.

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
On the membrane surface, there is a skin layer which have very fine thin pores and the objective materials would be separated according to its molecular size.

Hollow fiber shape is straw-like and thousands to several ten thousands of it is set inside the module.

Supplying the feed inside the hollow fiber is Inside-out filtration.

Suppling the feed outside the hollow fiber is Outside-in filtration.

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Widely prepared modules from lab. module for table test till large sized modules(F50type) for large scale facility.

品番 Model			分画分子量 Molecular weight cut-off	膜材質*1 Membrane material	膜面積 Membrane area (m <sup>2</sup> )	中空糸内径 Hollow fiber internal dia. (mm)	最高使用温度*2 Max. applicable temperature (°C)	pH範囲 pH range	初期ろ過水量 Initial pure water flux(Reference) (m <sup>3</sup> /Hr·0.1MPa)
モジュール型式 Module	膜種 Membrane	グレード Grade							
F50	FUC	1582	150,000	CA	50	0.8	40	4 ~ 8	28.5
F20	FUS	0382	30,000	PES	17	0.8	85	1 ~ 13	6.8
	FUC	1582	150,000		16		40	4 ~ 8	9.8
F10	FUY	03A1	30,000	PAN	4.2	1.0	45	2 ~ 11	1.2
	FUS	T653	6,000	PES	7.8	0.5	98	1 ~ 13	3.0
		0181	10,000		5	0.8			1.8
		0382	30,000		7.8	0.5			2.0
		0353			5	0.8			3.0
		1582	150,000		5	0.8			4.0
5082	500,000	5	0.8	3.9					
F03	FUS	1041	100,000	PES	2.3	0.4	98	1 ~ 13	1.0
	FUS	0353	30,000		2.1	0.5			0.7
		0382			1.4	0.8			0.6

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※1 CA = cellulose acetate, PES = polyethersulfone, PAN = polyacrylonitrile  
 ※2 It may differ depending on operating conditions. Especially, in case of high temperature use always please consult in advance.

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■ Most suitable for solid-liquid separation of high viscosity and turbidity  
Due to wider cut area of feed channel, it is unlikely to be fouled.

■ Effective for recovery of valuable matters and application of foods which are spoil easily.  
Simple structure of conc. water is unlikely to be retained.

■ The low elusion risk of module parts to the products.  
Because the glue is not used for module parts, so no influence on the flavor and smell through the process.

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

The tubular type membrane module is a cylindrical membrane module with a diameter of 10 to 15 mm and a length of around 3 m.

The filtration method is basically an internal pressure filtration method, and water or liquid that has passed through the membrane passes through the cylindrical support pipe and flows out through the hole of the pressure-resistant tube.

Compared to hollow fiber membrane modules and spiral membrane modules, the structure is simple and the raw water flow path is large, so it is characterized by being difficult to be fouling even if the suspended matter is directly filtered.

Also, with a membrane module structure suitable for concentration and separation, only cylindrical support tube can be replaced.

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Series type: MH, P18A, P18B, M4

Raw water flows through elements arranged in series type.

Parallel type: P18LP

Raw water flows in the elements arranged in the parallel type.

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■ Most suitable for concentration of saccharide, amino acid, peptide from ion or at manufacturing process. (RO membrane)

■ Most suitable for refinery of dye, water solute polymer. (NF membrane)

■ High maintainability, compact  
Compact size due to high density membrane structure makes membrane replacement work easier.

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Spiral type membrane module has a structure of flat sheet type separation membrane in an envelope shape. Membrane element would be set in the vessel (pressure resistant vessel) and used.

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
In membrane treatment, no heat is required to concentration and separation.  
It is very useful treatment method of processing the food sensitive for heat  
with advantage of removing the unwanted components without change of phase.

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
Refinery and concentration of functional food by membrane.

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Membrane module available with continuous operation for refinery and concentration have been widely introduced at the producing process of the dye, pigment, water soluble polymer. The needs of refinery and concentration by membrane process is growing because of recent development of the nanosized technology.

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- Demineralization and refinery of dyes ▪ Refinery and concentration of water-soluble polymers
- Refinery and concentration of metal nanoparticles ▪ Concentration of latex

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
Refinery, concentrate or recover large molecular weight substances

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※ Diafiltration method, diafiltration  
A method of purifying a target substance by diluting the feed solution with a solvent (mainly water) or diluting it beforehand and filtering the solvent and low molecular weight substance by filtration.  
(Quoted from JIS K 3802: 2015)

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With high separation performance and heat and chemical resisting character, widely introduced to the fields from refinery water/injection water production at pharmaceutical to cell separation/concentration of enzymes.

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
- Enzyme refinery and concentration ▪ Purified water production ▪ Injection water production
- Pyrogen in preparation process ▪ Protein
- Concentrate and desalted purification of enzyme preparation

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※1 Purified water is refined water by distillation, ion exchange, reverse osmosis, ultra filtration or combination of them.

※2 Water for Injection is produced by distillation or ultrafiltration of purified water. These process are defined by Japan Pharmacopeia

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Because water treatment by membrane filtration method can remove the bacteria and virus physically, the treated water is safe and high quality.

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
■ 従来法との比較 Comparison with the conventional method

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
- With the introduction of membrane treatment methods, the simplification of the drinking water plant can be possible. Thanks to this, shortening the building period of the drinking water plant facility and space saving can be also realized.

従来法(凝集沈殿~砂ろ過) Conventional method (coagulation sedimentation - sand filtration)

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The ground water can be safe and tasty water with ultrafiltration membrane (UF membrane) treatment. This treatment have been adopted at diverse facilities and companies.

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- Turbidity removal, disinfection
- Water and sewage water cost cut.

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- Improvement of drinking water safety, quality improvement of manufacturing water
- Improvement of product quality, water and sewage water cost cut.

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Turbidity removal and sterilized of sea water can be done with ultrafiltration membrane (UF membrane) treatment. This treatment have been adopted at marine products company, fishing port and farm. And also the sea water can be desalt with reverse osmosis (RO membrane) treatment.

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For cost reduction and recovering the values in the waste water, our membrane products have been widely adopted at the factories and offices. We have a lot of references at night soil treatment plants with our tubular type membrane module used as side stream MBR.

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- Separation of water-soluble oil
- Wafer polishing drain
- Reduction of waste water
- Night soil treatment
- Recycling through advanced treatment of waste water(grey water)

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We are the only manufacturer in Japan that manufactures membrane raw materials,  
manufacturing modules, processing systems, consistently.

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MOLSEP® is a trademark of Daisen Membrane-Systems' separation membrane module and system.  
Taking full advantage of the technology and abundant achievements of separation membrane modules at Daicel Co., Ltd.,  
we will provide total support from laboratory testing to engineering and maintenance.

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We will propose the combinations of membrane modules  
and pre-treatment.

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Covering wide water treatment fields including drinking water,  
food, chemical, waste water.

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We can offer the maintenance work after installation.

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We propose optimum membrane processing method according to customer's needs.

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You can send the samples to us and we will have a lab. test.  
Or also you can have the test by yourself.

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We will simulate initial cost and running cost.  
We will also propose processing flow and installation layout.

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We can have a running test at the site and then narrow the operational conditions.

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We can design and assemble the membrane treatment system.

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After installation of membrane treatment system, we will have  
a trial running and adjustment at the site.

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Our specialized staff will be responsible for the maintenance work of  
membrane treatment system including membrane replacement work,  
chemical replenishing work, routine check of the system, water quality analysis.